



## Pulse oximeter S250

### User manual



Version: 2.0  
Date: 2020-07-08



Thank you for choosing this product from Scaleson US (haftungsbeschränkt). It is equipped with all the features of the latest technology and has been optimised for easy operation.

Please read the instructions for use carefully before use.

If you have any questions or if problems occur with your device that are not covered in the instruction manual, please contact your Scaleson US (haftungsbeschränkt) service partner or our customer service department:

Phone: +49 3377 / 9979 395  
E-Mail: service@scaleson.net

This manual describes the commissioning and operation of the Scaleson Pulse Oximeter S250.

**Scope of delivery:**  
Pulse oximeter  
2 x AAA alkaline batteries  
Lanyard  
User manual

**Intended use:**  
The pulse oximeter is suitable for private use (at home) as well as for use in the medical sector (hospitals, medical establishments). Only use the pulse oximeter on humans to measure the arterial oxygen saturation (SpO<sub>2</sub>) of haemoglobin and the heart rate (pulse rate).

### 1 Warnings and safety notes

Read these instructions for use carefully. Non-observance of the following information may result in personal injury or material damage. Store these instructions for use and make them accessible to other users. Make sure you include these instructions for use when handing over the device to third parties.

Please report all serious incidents that have occurred in connection with the product directly to the manufacturer and the competent authority in your country.

**⚠ WARNING**

- Check to ensure that the package contains all the parts that should be included in the delivery.
- Check the pulse oximeter regularly before use to ensure that there is no visible damage to the device and the batteries are still sufficiently charged. In case of doubt, do not use the device and contact the customer services or an authorised retailer.
- Do not use any additional parts that are not recommended by the manufacturer or offered as equipment.
- Under no circumstances should you open or repair the device yourself, as faultless functionality could no longer be guaranteed thereafter. Failure to comply will result in voiding of the warranty. For repairs, please contact the customer services or an authorised retailer.
- Using the device for long periods may cause pain for people with circulatory disorders. Therefore do not use the pulse oximeter for longer than approx. 2 hours on one finger.
- The pulse oximeter displays a current measurement but cannot be used for continuous monitoring.
- The pulse oximeter does not have an alarm function and is therefore not suitable for evaluating medical results.

- Do not self-diagnose or self-medicate on the basis of the measurements without consulting your doctor. In particular, do not start taking any new medication or change the type and/or dosage of any existing medication without prior approval.
- Do not look directly inside the housing during the measurement. The red light and the invisible infra-red light in the pulse oximeter are harmful to your eyes.
- This device is not intended for use by people (including children) with restricted physical, sensory or mental skills or a lack of experience and/or a lack of knowledge, unless they are supervised by a person who has responsibility for their safety or they receive instructions from this person on how to use the device. Children should be supervised around the device to ensure they do not play with it.
- Neither of the displays for the pulse wave and pulse bar allows the strength of the pulse or circulation to be evaluated at the measurement site. Rather, they are exclusively used to display the current visual signal variation at the measurement site and do not enable reliable diagnostics for the pulse.

Do NOT use the pulse oximeter

- if you are allergic to rubber products.
- if the device or the finger you are using is damp.
- on small children or babies.
- during an MRI or CT scan.
- whilst taking a blood pressure measurement on the same arm using a cuff.
- on fingers that have nail varnish on, are dry or have a plaster or other dressing on them.
- on large fingers that do not fit into the device easily
- on fingers with anatomical changes, oedemas, scars or bumps.
- on fingers that are too small, as with small children for example.
- on patients who are not steady at the site of application (e.g. trembling).
- near flammable or explosive gas mixtures.

Non-observance of the following instructions can lead to incorrect or failed measurements.

- There must not be any nail varnish, artificial nails or other cosmetics on the finger to be measured.
- Ensure that the finger nail on the finger to be measured is short enough that the fingertip covers the sensor element in the housing.
- Keep your hand, finger and body steady during the measurement.
- For people with cardiac arrhythmia, the measurement values of SpO<sub>2</sub> and the heart rate may be incorrect or the measurement may not be possible at all.
- In cases of carbon monoxide poisoning, the pulse oximeter displays a measurement value that is too high.
- To avoid falsifying the measuring result, there should not be any strong light sources (e.g. fluorescent lamps or direct sunlight) in the immediate vicinity of the pulse oximeter.
- People with low blood pressure, who suffer from jaundice or take medication for vascular contraction, may experience incorrect or falsified measurements.
- Incorrect measurements are likely for patients who have been administered medical eye in the past or for those who have abnormal haemoglobin levels. This applies in particular for cases of carbon monoxide poisoning and methaemoglobin poisoning, which can occur for example from the administration of local anaesthetics or from an existing methaemoglobin reductase deficiency.
- Protect the pulse oximeter from dust, shocks, moisture, extreme temperatures and explosive materials.

### 2 Device description




Fig. 1: Front view

### 3 Technical data

Model: YM301

Classification: IP22, application part type BF

Display type: OLED display

Display colour: blue

Displayed values: SpO<sub>2</sub>, Pulse Rate, P1%, Bar Graph Battery Indicator, Pulse Wave

Measurement method: Non-invasive measurement of arterial oxygen saturation of haemoglobin and pulse rate in finger

Measurement range: SpO<sub>2</sub> 90 – 99 %  
Pulse 25 – 250 beats/minute

Accuracy: SpO<sub>2</sub> 70 – 100 %: ± 2 %  
SpO<sub>2</sub> 0 – 69 %: unspecified  
Pulse 0 – 250 bpm, ± 3 beats / minute

Dimensions: 57 x 30 x 31 mm (LxBlxH)

Weight: 28 g (batteries not included)

Power supply: Two 1.5 V AAA (LR03) batteries

Permissible operating conditions: +10 °C to +40 °C, 10 to 95% relative humidity (non-condensing), 700 – 1060 hPa ambient pressure

Permissible storage conditions: -20 °C to +60 °C, 0 to 95% relative humidity (non-condensing), 500 – 1074 hPa ambient pressure

**Note:**

1) The claim for oxygen saturation accuracy should be supported by clinical studies covering the entire claimed range. The fraction of inspired oxygen (FiO<sub>2</sub>) delivered to test subjects is varied to achieve a series of targeted steady-state saturation periods over the specified SpO<sub>2</sub> accuracy range (e.g. 70 % to 100 %), then the SpO<sub>2</sub> accuracy is calculated by comparing SpO<sub>2</sub> readings of the pulse oximeter to the values of SaO<sub>2</sub> determined with a Co-Oximeter.

2) The clinical trial included 11 subjects of different skin colour, including 6 males and 5 females, with an age range of 18 to 46 years.

### 4 Commissioning

#### 4.1 Inserting the batteries

Insert the two AAA batteries into the battery compartment with the correct polarity (see fig. 2)

Close the battery compartment lid again.




Fig. 2: Inserting the batteries

#### 4.2 Attaching the lanyard

Push the narrow end of the lanyard through the holder.

Tighten the other end of the lanyard through the loop of the narrow end. (see Fig. 3)




Fig. 3: Attaching the lanyard

### 5 Operation

Insert one finger into the finger opening of the pulse oximeter as shown (see Fig. 4) and hold it steady.

Press the power button. The pulse oximeter begins its measurement. Do not move during the measurement.

Your measurement values will appear on the screen after a few seconds.




Fig. 4: Optimal finger positioning

**Note**

When you remove your finger from the pulse oximeter, the device will automatically switch off after approx. eight seconds.

**Function button**

The function button on the pulse oximeter has two functions in total:

- Switch on function: When the pulse oximeter is switched off you can hold down the function button briefly to switch it on.
- Setting the parameters: To set the individual parameters, such as the display brightness, press and hold the function key during operation.

### Setting the parameters

After switching on the oximeter, press and hold the function key for about 2 seconds. The oximeter calls up the interface for setting the parameters.




Fig. 5: Interface 1 and 2 settings menu

Press the function key briefly to cycle through the individual parameters. The star shows which parameter you are currently at. By pressing and holding the function key you can select the parameter and then make the desired changes. To exit the setting menu, select the item "Exit" and keep the function key pressed for a long time.

Display	Parameter	Setting options
Demo	Demo view	By activating the demo view, a sample view of the measured values appears permanently on the display.
Restore	Factory setting	By pressing the function key for a long time, you can reset all settings to factory default.
Brightness	Display brightness	Press the button and select the brightness level in the range from 0 to 5. The larger the value, the brighter the screen.
SpO2 Alm hi	Upper limit of the SpO2 alarm	Press the button and select the upper limit of the SpO2 alarm: 95, 96, 97, 98, 99 or 100.
SpO2 Alm low	Lower limit of the SpO2 alarm	Press the button and select the lower limit of the SpO2 alarm: 80, 81, 82, 83, 84, 85, 86, 87, 88, 89, 90, 91, 92, 93 or 94.
Pr Alm hi	Upper limit of the pulse alarm	Press the button and select the upper limit of the pulse alarm: 90-200 beats/minute.
Pr Alm low	Lower limit of the pulse alarm	Press the button and select the lower limit of the pulse alarm: 40-200 beats/minute.
+/-	Increase or decrease the limit values	Spikes "+" to increase the limit values listed above. Select "-" to decrease the limit values.

### 6 Maintenance and cleaning

**⚠ IMPORTANT:**  
Do not use high-pressure sterilisation on the pulse oximeter.

The recommended disinfectant is 70% ethanol. Never use ETO or formaldehyde for disinfection.

Under no circumstances should you hold the pulse oximeter under water, as this can cause liquid to enter and damage the pulse oximeter.

- Clean the housing and the interior rubber surface with a soft cloth dampened with medical alcohol after each use.
- If a low battery status appears on the display of the pulse oximeter, change the batteries.
- If you are not going to use the pulse oximeter for more than one month, remove both batteries from the device to avoid possible leaking.

### 8 Disposal

For environmental reasons, do not dispose of the device in the household waste at the end of its useful life. Dispose of the unit at a suitable local collection or recycling point.

Dispose of the device in accordance with EC Directive – WEEE (Waste Electrical and Electronic Equipment).

If you have any questions, please contact the local authorities responsible for waste disposal.

The empty, completely flat batteries must be disposed of through specially designated collection boxes, recycling points or electronics retailers. You are legally required to dispose of the batteries.

The codes below are printed on batteries containing harmful substances:

Pb = Battery contains lead,  
Cd = Battery contains cadmium,  
Hg = Battery contains mercury.

### 9 Malfunctions - causes and elimination

Malfunction	Possible causes	Description / solution
The SpO <sub>2</sub> and PR values are not displayed normally	The measuring finger is too big or too small. Excessive ambient light. User's blood perfusion is very low.	Select a suitable finger. Avoid excessive exposure to ambient light. Warm the finger and try again.
The display is suddenly off	The unit switches off automatically after 8 sec. if it does not receive any signals.	Normal state. No error.
The SpO <sub>2</sub> and PR values are not displayed stably	The batteries are discharged or almost discharged.	Replace batteries.
	Finger, hand or body is in motion.	Keep fingers, hand and body still during measurement.
	Not used in the working environment recommended in this manual.	Please use in a normal working environment.
	The device is damaged.	Please contact the dealer.

### 6 Note on electromagnetic compatibility

The device complies with the requirement of EN 60601-1-2 "Electromagnetic Compatibility - Medical Electrical Equipment".

Immunity test	IEC 60601-1 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD)	± 8 kV contact discharge ± 2, ± 4, ± 8, ± 15 kV air discharge	± 8 kV contact discharge ± 2, ± 4, ± 8, ± 15 kV air discharge	Floors should be made of wood or concrete or covered with electrostatic conductive material, the relative humidity must be at least 30%.

Electrical fast transient/burst emissions	N/A	N/A	N/A
7 IEC 10100-3-3			


Surge	± 1.5 kV, ± 1 kV differential mode line-line	N/A	N/A
IEC 61000-4-5			

Power frequency (50/60 Hz) magnetic field	30 A/m, 50/60 Hz	30 A/m, 50/60 Hz	Magnetic fields at the mains frequency should correspond to the typical values found in the business and hospital environment.
IEC 61000-4-4			

NOTE: UT is the AC mains voltage before applying the test level.

Guidelines and manufacturer's declaration - electromagnetic immunity	IEC 60601-1 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-3	3 Vrms, 150 kHz to 80 MHz 6 Vrms, 150 kHz to 80 MHz outside ISM bands	N/A	Portable and mobile RF communication equipment should be used at no less distance from the pulse oximeter, including the hands, than the recommended separation distance, calculated according to the separation frequency.
Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range.			
c) The field strength of stationary transmitters, such as base stations of radio telephones and mobile land mobile, amateur radio stations, AM and FM radio and television transmitters cannot be predicted theoretically with any accuracy. To determine the electromagnetic environment with regard to stationary transmitters, a study of the site should be considered.			

Interference may occur in the vicinity of equipment marked with the following symbol:




NOTE 1: At 80 MHz and 800 MHz the higher frequency range applies. These guidelines may not be applicable in all cases. The propagation of electromagnetic quantities is influenced by absorptions and reflections of buildings, objects and people.

NOTE 2: If the measured field strength at the site where the equipment is used exceeds the above compliance level, the equipment should be observed to verify that it functions as intended. If unusual performance characteristics are observed, additional measures may be necessary, such as changing or relocating the equipment.

d) Over the frequency range from 150 kHz to 80 MHz, the field strength should be less than 3 V/m.

If the measured field strength at the site where the equipment is used exceeds the above compliance level, the equipment should be observed to verify that it functions as intended. If unusual performance characteristics are observed, additional measures may be necessary, such as changing or relocating the equipment.

d) Over the frequency range from 150 kHz to 80 MHz, the field strength should be less than 3 V/m.



Scaleson US (haftungsbeschränkt)  
Mellnerstr. 21  
15806 Zossen  
Telefon +49 3377 / 9979 395  
info@scaleson.net  
www.scaleson-shop.net

Shenzhen Yim Lik Technology Co., Ltd  
Address: 205, Building A,  
Tinghe Industrial Park,  
Changshengqian Street, Longde Village,  
Pinghuang District, 518118 Shenzhen  
Made in China.

Shanghai International Holding Corp. GmbH (Group)  
Address: Effelrodweg 86, 20527 Hamburg, Germany  
E-Mail: a.holding@international.com

© Scaleson US (limited liability), publication, reproduction and any form of use, as well as disclosure to third parties – even in part or in revised form – without the consent of Scaleson US (limited liability) is prohibited. Technical changes reserved.  
470 080 012 (Version 1.0) Stand 08/2025