

BLOOD PRESSURE MONITOR S260.10



Brand: Scaleson

- Precise measurement in accordance with WHO standards
- For 2 users – 120 memory slots each
- Large LED display
- Voice guidance in German and English
- Smart pressure regulator & comfortable measurement
- USB Type-C port
- App functionality



Product Description

Our fully automatic digital blood pressure monitor offers you an accurate and convenient way to monitor your blood pressure regularly – both at home and in healthcare settings. Using advanced oscillometric technology, the device reliably measures systolic and diastolic blood pressure as well as your pulse.

Measurements are taken conveniently on the upper arm. The ergonomically designed cuff is suitable for a circumference of **22 to 40 cm** and ensures comfortable use for adults and young people aged 12 and over.

High-performance technology for reliable results

- Fully automatic measurement of **blood pressure and pulse**
- Precise results thanks to **oscillometric technology**
- **Six-segment scale display** for clear classification of readings
- Easy-to-read LED display
- Easy to use thanks to intuitive button controls

The device impresses with its well-thought-out design: it features a compact housing, a cuff connection, a connection for the AC adapter and a clear, user-friendly display. The built-in rechargeable battery allows for flexible use – at home or in a professional setting.

Easy to use

For best results, we recommend fully charging the device before using it for the first time. It will then be ready for use at any time and will reliably assist you in regularly monitoring your blood pressure.

Important note

This blood pressure monitor is intended for measuring blood pressure in people aged 12 and over. The results should only be interpreted by qualified medical professionals.

Rely on modern measurement technology and a well-thought-out design – for greater safety and control in everyday life.

Wide range of applications...



At home



Pharmacy



Doctor

Perfect for reliable measurements in the comfort of your home, this blood pressure monitor offers consistent accuracy. Two users can operate the device independently without any issues. Individual user profiles ensure that results are clearly organized and easy to distinguish.

Designed specifically for pharmacy use, this device supports a wide range of everyday applications. Whether measuring patients' blood pressure, offering advice, or demonstrating its use – it performs efficiently in every situation. Quick and precise results are ensured thanks to its simple, automatic operation. Additionally, the built-in voice function makes the process easier and improves overall user-friendliness, especially for older users.

In doctors' practices, this blood pressure monitor proves highly reliable for daily use. Fast and accurate readings support efficient patient assessments. Healthcare professionals benefit from its dependable performance during routine care.

The special feature of our forehead and ear thermometer S240...



Voice function

The built-in voice function provides clear guidance through the measurement process and, if desired, reads the results aloud – a particularly convenient feature for older people. Via the app connection, all measurement data can be conveniently saved, analysed and tracked over time on your smartphone



WHO scale

The built-in WHO scale shows at a glance how the measured blood pressure should be interpreted. The colour-coded display is based on World Health Organisation (WHO) guidelines and indicates whether the blood pressure is normal, elevated or too high. This makes it easy to understand the readings straight away



App feature

This blood pressure monitor is ideal for use in doctors' surgeries. It reliably delivers fast and accurate readings in day-to-day practice. The device also efficiently supports healthcare professionals in their daily patient care.

[Download User instructions](#)

[Download Quick](#)

Technical Data

Order no.	S260.10.000
Measuring range	Taille : très grande Poids de l'article : 381 grammes Dimensions de l'article : 17 L x 8,5 l x 12,6 P cm
Memory capacity	2 x 120 emplacements mémoire
Power supply	Batterie rechargeable et USB-C
Batteries	Alimentation : batterie rechargeable

Weight	381 grammes
Dimensions	17 cm (L) x 8,5 cm (l) x 12,6 cm (P)