


# ANIMAL SCALE S230.15



Brand: Scaleson

- 800 x 500 mm
- Small animal scales
- Removable rubber mat
- Maximum load 100 kg, 10 g increments
- Stable base with rubber feet



# Product Description

The digital animal scale impresses with its ease of use, high accuracy and versatile functions. With a load capacity of up to 100 kg and a resolution of 10 g, it is ideal for everyday weighing of animals.

The included rubber mat protects the surface of the scale and ensures that dogs, cats, rabbits and many more can be weighed without fear.

All animals weighing up to 100 kg are welcome on these scales.  
By the way: small movements do not affect the measurement results.

Kilograms (kg), pounds (lb) or ounces (oz) – select the unit that suits you best at the touch of a button.

## Wide range of applications...



### Veterinary hospital

The S230 animal scale makes weighing animals easy, safe and reliable. It allows weight data to be recorded quickly and stress-free. The non-slip surface ensures that animals feel comfortable and remain on the scale without fear. Even slight movements do not affect the measurement, so reliable results are available at all times. The appropriate weight unit can be set at the touch of a button. The scale combines precision, comfort and flexibility – a valuable aid for vets and animal owners.



### Animal shelter

The S230 animal scale makes weighing animals easy, safe and reliable. It allows weight data to be recorded quickly and stress-free. The non-slip surface ensures that animals feel comfortable and remain on the scale without fear. Even slight movements do not affect the measurement, so reliable results are available at all times. The appropriate weight unit can be set at the touch of a button. The scale combines precision, comfort and flexibility – a valuable aid for vets and animal owners.



### At home

The S230 animal scale makes daily work easier for animal shelters by weighing animals quickly and reliably. This means that staff always have reliable data at their fingertips to keep track of the health of their charges without much effort. Whether for admitting new animals or checking their nutritional status, the S230 ensures efficiency and safety in the workflow.


## What makes our S230 animal scales special...



Easy to use



Non-slip rubber mat



Stable footing



Use as a small animal scale

The scale can be switched on and off very easily using the 'ON/TARE/OFF' button. The tare function allows you to determine the net weight, e.g. without blankets or towels underneath. The unit button 'UNIT' can be used to change the unit from kg to lbs. By pressing the 'MEMORY' button, you can display the last weight measured at any time.

The highlight of the Scaleson S230 small animal scale is the removable rubber mat. The animals do not slip and have a firm grip during the weighing process.

The underside of the scales is equipped with four rubber feet. This ensures a secure hold on a wide variety of surfaces.

The Scaleson S230 small animal scale is suitable for all types of small animals and impresses with its ease of use, high accuracy and versatile functions. Small movements do not affect the measurement result, enabling accurate weight determination at all times.

[Download Manual](#)

## Technical Data

|                      |                                                     |
|----------------------|-----------------------------------------------------|
| <b>Artikelnummer</b> | S230.05.100                                         |
| <b>Dimensions</b>    | 800 x 500 x 40 mm                                   |
| <b>Maximum load</b>  | 100 kg                                              |
| <b>Division</b>      | 10 g                                                |
| <b>Batterien</b>     | 4 x AAA-Batterien (nicht im Lieferumfang enthalten) |