

BLOOD PRESSURE MONITOR S260.20



Brand: Scaleson

- Accurate measurement in accordance with WHO standards
- Large LED display
- Voice guidance in German and English
- Smart pressure regulator & comfortable measurement
- USB Type-C port
- App functionality
- Includes a handy storage compartment for the cuff



Product Description

Keep an eye on your health – quickly, easily and accurately. Our fully automatic digital blood pressure monitor allows you to conveniently measure your blood pressure and pulse right in the comfort of your own home.

Thanks to state-of-the-art measurement technology, the device reliably measures your systolic and diastolic blood pressure as well as your pulse – simply at the touch of a button. The soft upper-arm cuff (22–40 cm) ensures a comfortable fit and straightforward use.

The easy-to-read display clearly shows your readings. An integrated colour scale based on WHO standards helps you to better interpret and quickly understand your results.

Your benefits:

- Simple, automatic measurement at the touch of a button
- Reliable results for blood pressure and pulse
- Intuitive display with colour scale for guidance
- Comfortable upper arm cuff (22–40 cm)
- Ideal for regular monitoring at home

This blood pressure monitor is suitable for users aged 12 and over and helps you keep a regular eye on your readings – for greater peace of mind in everyday life.

Wide range of applications...




At home

This device is the ideal choice for reliable blood pressure monitoring in the comfort of your own home. Its simple design makes it easy for two people to use. Thanks to customisable user profiles, all readings are clearly organised and kept separate from one another.



Pharmacy

Especially in the current situation, it is becoming more and more important to check regularly in the office or at the workplace whether employees or external visitors may have a fever. With the Scaleson forehead and ear thermometer S240 you can check this contactlessly in just a few seconds.



Doctor

In everyday clinical practice, the device impresses with its reliable and rapid measurement performance. It supports healthcare professionals in their daily patient care and provides accurate blood pressure readings in a short space of time. This ensures an efficient workflow.

The special feature of our forehead and ear thermometer S260.20...



Voice function

The built-in voice guidance provides clear instructions throughout the measurement process and, if desired, reads the results aloud – a particularly convenient feature for older people. Via the app connection, all measurement data can be conveniently saved, analysed and tracked over time on your smartphone.



WHO scale

The built-in WHO scale shows at a glance how the measured blood pressure should be classified. The colour-coded display is based on World Health Organisation (WHO) guidelines and indicates whether the blood pressure is normal, elevated or high. This makes it easy to understand the readings straight away.



App feature

Using the right app, you can easily save and record all your blood pressure readings. This allows you to clearly track your readings over time, analyse them in the long term and, if necessary, easily share them with your doctor or pharmacist.

[Download User instructions](#)

[Download Quick](#)

Technical Data

Order no.	RF-1T6Z-RN6Q
Measuring range	Size: Extra large Item weight: 365 grams Item dimensions: 13L x 10W x 16H cm
Memory capacity	2x120 memory slots
Power supply	Battery-powered, requires 3 x AA batteries (not included), USB Type-C port
Batteries	3 x AA batteries (not included)
Weight	365 grams
Dimensions	163,3 x 39,2 x 38,9 mm